## **JAVALYTICS**<sup>TM</sup>

by MADISON INSTRUMENTS, INC.

## JAV-RDA-H Handheld Degree of Roast Analyzer



## Portable, battery powered, simple to operate, and affordable!

## **Key Specifications:**

#### **Power:**

Lithium polymer battery operated, charged in a convenient charging cradle. Fully charges in approximately 2 ½ hours for about 2 hours of continuous operation. The AC adapter has 5VDC, 2A output.

## Weight:

Analyzer: 5.5 ounces (156 grams)	Charging cradle: 3.5 ounces (98 grams)
AC adapter: 3.3 ounces (93 grams)	Cal Disk: 0.7 ounces (21 grams)
Sample Tray: 3.4 ounces (95 grams)	Measuring cup: 0.4 ounces (12 grams)

**History:** 100 records stored.

**Temperature:** Storage: -10 - 60C. Operation: 0 - 50C.

**Compliance:** CE and part 15 of FCC rules.


4319 Twin Valley Rd Ste 1, Middleton, WI 53562-4415 Toll Free in the USA: 888-349-3847 Phone: +1 608-833-1854 Fax: +1 608-833-1859

EMAIL: paulh@javalytics.com WEB: www.javalytics.com

## **JAVALYTICS**<sup>TM</sup>

by MADISON INSTRUMENTS, INC.

## Booting up



## Main Menu



Choose calibration (CAL), measurement (MSR), review historical records (REC), or setup (SET).

### Measurement



Choose Ground or Whole bean. Whole bean will prompt through three samples for averaging.

## Then...



# Position and click OK.



#### And the result:



4319 Twin Valley Rd Ste 1, Middleton, WI 53562-4415 Toll Free in the USA: 888-349-3847 Phone: +1 608-833-1854 Fax: +1 608-833-1859